


SURGICRYL[®]

MONOFILAMENT

POLYDIOXANONE



Type	Monofilament
Composition	Polydioxanone, a polymer made from polyester poly (p-dioxanone)
Coating	None
Colour	Violet
Tissue reaction	Minimal
Absorption	The hydrolytic action by which the material is broken down results in total absorption in approximately 180 to 210 days. Approximately 50% of tensile strength remain after 42 days.
Presentation	Box with 12 sutures



*For quality
&
safety*



www.sutures.be

SURGICRYL[®]

MONOFILAMENT

POLYDIOXANONE

CHARACTERISTICS

- High and long-lasting tensile strength
- Synthetic
- Breaks down by hydrolysis
- Monofilament with microscopic uniform diameter
- Pliability and softness
- Hermetically sealed packing

BENEFITS

- Outstanding strength retention for extended wound support
- Minimal tissue reaction
- Predictable absorption rate
- Soft passage through the tissues without sawing, tissue drag and trauma
- Absence of capillarity (no support of bacterial growth)
- Excellent knot tie down
- Good handling properties
- Easy and secure knotting
- Guaranteed seal and product sterility

INDICATIONS General soft tissue closing especially in


- General surgery
- Gastrointestinal surgery
- Ophthalmic surgery
- Orthopaedics
- Abdominal wall closure
- Gynaecology
- Plastic surgery
- Urology

CONTRAINDICATIONS Cardiovascular surgery
Neurosurgery
Prosthetic devices


Being absorbable not to use when wound support for more than 6 weeks is required.


Caution in patients with delayed wound healing (e.g. compromised immune system).

TENSILE STRENGTH



ABSORPTION



 Without needle	Metric	USP	Length	Order number
	5	2	150 cm	<i>violet</i> 13500015
	6	3&4	150 cm	<i>violet</i> 13600015

 Without needle	Metric	USP	Length	Order number
--------------------------------------------------------------------------------------------------	--------	-----	--------	--------------


Ligapack


















2 3/0 250 cm *violet* 13200025






 Round bodied taper point	Metric	USP	Length	Order number
------------------------------------------------------------------------------------------------------------	--------	-----	--------	--------------







1/2 circle

















HR 9 mm		0.7	6/0	75 cm	<i>violet</i> 13070109
		0.5	7/0	45 cm	<i>violet</i> 13050113
HR 13 mm		0.7	6/0	75 cm	<i>violet</i> 13070113
		1	5/0	75 cm	<i>violet</i> 13100113
		1.5	4/0	75 cm	<i>violet</i> 13150113
		2	3/0	75 cm	<i>violet</i> 13200113
HR 17 mm		0.7	6/0	75 cm	<i>violet</i> 13070117
		1	5/0	75 cm	<i>violet</i> 13100117
		1.5	4/0	75 cm	<i>violet</i> 13150117
HR 22 mm		2	3/0	75 cm	<i>violet</i> 13200117
		1.5	4/0	75 cm	<i>violet</i> 13150122
HR 26 mm		2	3/0	75 cm	<i>violet</i> 13200122
		3.5	0	75 cm	<i>violet</i> 13350126
HR 30 mm		1.5	4/0	75 cm	<i>violet</i> 13150126
		2	3/0	75 cm	<i>violet</i> 13200126
		3	2/0	75 cm	<i>violet</i> 13300126
HR 36 mm		3.5	0	75 cm	<i>violet</i> 13350130
		2	3/0	75 cm	<i>violet</i> 13200130
HR 40 mm		3	2/0	75 cm	<i>violet</i> 13300130
		4	1	90 cm	<i>violet</i> 13400130
		3.5	0	90 cm	<i>violet</i> 13350136
HR 48 mm		4	1	90 cm	<i>violet</i> 13400136
		5	2	90 cm	<i>violet</i> 13500136
HR 30 mm		3	2/0	75 cm	<i>violet</i> 13300140
		4	1	90 cm	<i>violet</i> 13400140
		5	2	90 cm	<i>violet</i> 13500140
HR 48 mm		3.5	0	90 cm	<i>violet</i> 13350148
		4	1	90 cm	<i>violet</i> 13400148
		5	2	90 cm	<i>violet</i> 13500148



 Round bodied taper point			Metric	USP	Length		Order number
3/8 circle 	DR 13 mm		1	5/0	75 cm	<i>violet</i>	13100213
			1.5	4/0	75 cm	<i>violet</i>	13150213
	DR 17 mm		1	5/0	75 cm	<i>violet</i>	13100217
			1.5	4/0	75 cm	<i>violet</i>	13150217



 Round bodied taper point double armed			Metric	USP	Length		Order number
1/2 circle 	2 x HR 13 mm		0.7	6/0	75 cm	<i>violet</i>	13070113-2
			1	5/0	75 cm	<i>violet</i>	13100113-2
	2 x HR 17 mm		1	5/0	90 cm	<i>violet</i>	13100117-2
			1.5	4/0	75 cm	<i>violet</i>	13150117-2
			2	3/0	75 cm	<i>violet</i>	13200117-2
	2 x HR 26 mm		2	3/0	90 cm	<i>violet</i>	13200126-2



 Round bodied taper point double armed			Metric	USP	Length		Order number
3/8 circle 	2 x DR 9 mm		0.7	6/0	75 cm	<i>violet</i>	13070209-2
	2 x DR 10 mm		0.5	7/0	75 cm	<i>violet</i>	13050210-2
	2 x DR 13 mm		0.5	7/0	75 cm	<i>violet</i>	13050213-2
			0.7	6/0	75 cm	<i>violet</i>	13070213-2
			1	5/0	90 cm	<i>violet</i>	13100213-2
	2 x DR 17 mm		0.5	7/0	75 cm	<i>violet</i>	13050217-2
			1.5	4/0	90 cm	<i>violet</i>	13150217-2



 Round bodied taper point loop	Metric	USP	Length	Order number	
1/2 circle 	 HR 40 mm	3.5	0	150 cm	<i>violet</i> 13350940
		4	1	150 cm	<i>violet</i> 13400940
	 HR 50 mm	3.5	0	150 cm	<i>violet</i> 13350950
		4	1	150 cm	<i>violet</i> 13400950
		5	2	150 cm	<i>violet</i> 13500950
	 HR 65 mm	4	1	150 cm	<i>violet</i> 13400965

 Reverse cutting	Metric	USP	Length	Order number	
3/8 circle 	 DS 12 mm	0.7	6/0	45 cm	<i>violet</i> 13071512
		1	5/0	45 cm	<i>violet</i> 13101512
		1.5	4/0	75 cm	<i>violet</i> 13151512
	 DS 16 mm	0.7	6/0	45 cm	<i>violet</i> 13071516
		1	5/0	75 cm	<i>violet</i> 13101516
		1.5	4/0	75 cm	<i>violet</i> 13151516
	 DS 16 mm Easy Pass	1	5/0	50 cm	<i>violet</i> 13101516P
	 DS 16 mm black	1	5/0	75 cm	<i>violet</i> 13101516VB
	 DS 19 mm	1	5/0	45 cm	<i>violet</i> 13101519
		1	5/0	75 cm	<i>violet</i> 13101519-75
		1.5	4/0	75 cm	<i>violet</i> 13151519
		1.5	4/0	90 cm	<i>violet</i> 13151519-90
		2	3/0	75 cm	<i>violet</i> 13201519
	 DS 19 mm black	1.5	4/0	75 cm	<i>violet</i> 13151519VB
	 DS 24 mm	1.5	4/0	75 cm	<i>violet</i> 13151524
2		3/0	75 cm	<i>violet</i> 13201524	
2		3/0	90 cm	<i>violet</i> 13201524-90	
3		2/0	75 cm	<i>violet</i> 13301524	
3.5		0	75 cm	<i>violet</i> 13351524	

 Reverse cutting		Metric	USP	Length		Order number
	DS 26 mm	2	3/0	75 cm	<i>violet</i>	13201526
		3	2/0	75 cm	<i>violet</i>	13301526
	DS 30 mm	3	2/0	75 cm	<i>violet</i>	13301530
		3.5	0	75 cm	<i>violet</i>	13351530
	DS 40 mm	4	1	90 cm	<i>violet</i>	13401540
		5	2	90 cm	<i>violet</i>	13501540
	DS 48 mm	4	1	90 cm	<i>violet</i>	13401548

 Reverse cutting		Metric	USP	Length		Order number
	HS 36 mm	3	2/0	75 cm	<i>violet</i>	13301436
		3.5	0	75 cm	<i>violet</i>	13351436
		4	1	75 cm	<i>violet</i>	13401436

 Round bodied taper cutting point		Metric	USP	Length		Order number
	HRT 26 mm	2	3/0	75 cm	<i>violet</i>	13200526
		3	2/0	75 cm	<i>violet</i>	13300526
	HRT 40 mm	3.5	0	90 cm	<i>violet</i>	13350540
		4	1	90 cm	<i>violet</i>	13400540

 Blunt loop		Metric	USP	Length		Order number
	HRB 50 mm	4	1	150 cm	<i>violet</i>	13401850