


SURGICRYL[®]

MONOFAST

POLYGLECAPRONE 25



Type	Monofilament
Composition	Polyglecaprone 25, a copolymer made of 75% glycolide and 25% ε-caprolactone
Coating	None
Colour	Violet / Beige - undyed
Tissue reaction	Minimal
Absorption	The hydrolytic action by which the material is broken down results in total absorption in approximately 90 to 110 days. Approximately 75% of tensile strength remain after 7 days.
Presentation	Box with 12 sutures



*For quality
&
safety*

The logo for smi, consisting of the lowercase letters "smi" in a bold, blue, sans-serif font with a white outline.

www.sutures.be

SURGICRYL[®]

MONOFAST

POLYGLECAPRONE 25

CHARACTERISTICS

- High tensile strength
- Synthetic
- Breaks down by hydrolysis
- Monofilament with microscopic uniform diameter

- Pliability and softness

- Hermetically sealed packing

BENEFITS

- Excellent tissue support
- Minimal tissue reaction
- Predictable absorption rate
- Soft passage through the tissues without sawing, tissue drag and trauma
- Absence of capillarity (no support of bacterial growth)
- Excellent knot tie down

- Good handling properties
- Easy and secure knotting

- Guaranteed seal and product sterility


INDICATIONS Superficial soft tissue and/or ligation
 especially in General surgery
 Skin closure
 Plastic surgery
 Gastrointestinal surgery
 Gynaecology
 Urology

CONTRAINDICATIONS Cardiovascular surgery
 Neurosurgery
 Ophthalmic surgery


Being absorbable not to use when an extended wound support is required.


Caution in patients with delayed wound healing (e.g. compromised immune system).


TENSILE STRENGTH





















ABSORPTION




















 Without needle	Metric	USP	Length		Order number
	3	2/0	150 cm	<i>violet</i>	16300015
	3.5	0	150 cm	<i>violet</i>	16350015
	4	1	150 cm	<i>violet</i>	16400015
	5	2	120 cm	<i>violet</i>	16500012
	5	2	150 cm	<i>violet</i>	16500015




 Without needle	Metric	USP	Length		Order number
Ligapack					
	5	2	250 cm	<i>violet</i>	16500025


 Round bodied taper point	Metric	USP	Length		Order number
HR 13 mm 	0.7	6/0	45 cm	<i>violet</i>	16070113
	1	5/0	45 cm	<i>violet</i>	16100113
HR 17 mm 	1.5	4/0	75 cm	<i>violet</i>	16150113
	1	5/0	45 cm	<i>violet</i>	16100117
HR 17 mm 	1.5	4/0	75 cm	<i>violet</i>	16150117
	2	3/0	75 cm	<i>violet</i>	16200117
HR 22 mm 	1.5	4/0	75 cm	<i>violet</i>	16150122
	2	3/0	75 cm	<i>violet</i>	16200122
HR 22 mm 	3	2/0	75 cm	<i>violet</i>	16300122
	1.5	4/0	75 cm	<i>violet</i>	16150126
1/2 circle  HR 26 mm 	2	3/0	75 cm	<i>violet</i>	16200126
	3	2/0	75 cm	<i>violet</i>	16300126
HR 26 mm 	3.5	0	75 cm	<i>violet</i>	16350126
	2	3/0	75 cm	<i>violet</i>	16200130
HR 30 mm 	3	2/0	75 cm	<i>violet</i>	16300130
	3.5	0	75 cm	<i>violet</i>	16350130
HR 30 mm 	4	1	75 cm	<i>violet</i>	16400130
	3	2/0	75 cm	<i>violet</i>	16300136
HR 36 mm 	3.5	0	90 cm	<i>violet</i>	16350136
	4	1	90 cm	<i>violet</i>	16400136
HR 40 mm 	3	2/0	75 cm	<i>violet</i>	16300140
	3.5	0	90 cm	<i>violet</i>	16350140
HR 40 mm 	4	1	90 cm	<i>violet</i>	16400140



 Round bodied taper point	Metric	USP	Length	Order number	
1/2 circle 		3.5	0	90 cm	<i>violet</i> 16350148
		4	1	90 cm	<i>violet</i> 16400148
		5	2	90 cm	<i>violet</i> 16500148
HR 65 mm		3	2/0	90 cm	<i>violet</i> 16300165





 Reverse cutting	Metric	USP	Length	Order number		
3/8 circle 		0.7	6/0	45 cm	<i>undyed</i> 16071512	
		0.7	6/0	75 cm	<i>undyed</i> 16071512-75	
		0.7	6/0	75 cm	<i>violet</i> 16071512V	
		1	5/0	45 cm	<i>undyed</i> 16101512	
		1	5/0	45 cm	<i>violet</i> 16101512V	
		1.5	4/0	45 cm	<i>undyed</i> 16151512	
	DS 12 mm Easy Pass		0.7	6/0	50 cm	<i>undyed</i> 16071512P
			1	5/0	50 cm	<i>undyed</i> 16101512P
			1.5	4/0	50 cm	<i>undyed</i> 16151512P
	DS 16 mm		1	5/0	75 cm	<i>undyed</i> 16101516
			1	5/0	75 cm	<i>violet</i> 16101516V
			1.5	4/0	75 cm	<i>undyed</i> 16151516
			1.5	4/0	75 cm	<i>violet</i> 16151516V
			2	3/0	75 cm	<i>undyed</i> 16201516
	DS 16 mm Easy Pass		1	5/0	50 cm	<i>undyed</i> 16101516P
			1.5	4/0	50 cm	<i>undyed</i> 16151516P
	DS 16 mm black		1	5/0	75 cm	<i>violet</i> 16101516VB
	DS 19 mm		1	5/0	45 cm	<i>undyed</i> 16101519
			1	5/0	75 cm	<i>violet</i> 16101519V
			1.5	4/0	75 cm	<i>undyed</i> 16151519
1.5			4/0	75 cm	<i>violet</i> 16151519V	
2			3/0	75 cm	<i>undyed</i> 16201519	
2			3/0	75 cm	<i>violet</i> 16201519V	
3			2/0	75 cm	<i>undyed</i> 16301519	
DS 19 mm Easy Pass		1.5	4/0	50 cm	<i>undyed</i> 16151519P	
		2	3/0	50 cm	<i>undyed</i> 16201519P	

Reverse cutting		Metric	USP	Length		Order number
3/8 circle 	DS 19 mm black 	1	5/0	75 cm	<i>violet</i>	16101519VB
		1.5	4/0	75 cm	<i>violet</i>	16151519VB
	DS 24 mm 	2	3/0	75 cm	<i>undyed</i>	16201524
		3	2/0	75 cm	<i>undyed</i>	16301524
	DS 24 mm Easy Pass 	2	3/0	75 cm	<i>undyed</i>	16201524P-75
		3	2/0	75 cm	<i>undyed</i>	16301524P-75
	DS 26 mm 	1.5	4/0	75 cm	<i>undyed</i>	16151526
		2	3/0	75 cm	<i>undyed</i>	16201526
		2	3/0	75 cm	<i>violet</i>	16201526V
		3	2/0	75 cm	<i>undyed</i>	16301526
		3	2/0	75 cm	<i>violet</i>	16301526V
	DS 26 mm Easy Pass 	3	2/0	75 cm	<i>undyed</i>	16301526P-75
	DS 30 mm 	3	2/0	75 cm	<i>undyed</i>	16301530
		3.5	0	75 cm	<i>violet</i>	16351530V
DS 40 mm 	3.5	0	75 cm	<i>undyed</i>	16351540	
	4	1	90 cm	<i>undyed</i>	16401540	

Reverse cutting		Metric	USP	Length		Order number
1/2 circle 	HS 36 mm 	3	2/0	75 cm	<i>undyed</i>	16301436
		3	2/0	75 cm	<i>undyed</i>	16301436V
		3.5	0	75 cm	<i>undyed</i>	16351436
	HS 40 mm 	4	1	90 cm	<i>undyed</i>	16401440

Cutting		Metric	USP	Length		Order number
Straight 	GS 60 mm	2	3/0	75 cm	<i>undyed</i>	16202160
		3	2/0	75 cm	<i>undyed</i>	16302160

Round bodied taper cutting point		Metric	USP	Length		Order number
1/2 circle 	HRT 40 mm 	3.5	0	90 cm	<i>undyed</i>	16350540

 Blunt		Metric	USP	Length		Order number
1/2 circle 	HRB 40 mm		3.5	0	90 cm	<i>violet</i> 16351840
	HRB 48 mm		4	1	90 cm	<i>violet</i> 16401848