

Consumable Catalogue

N2 (UK) Ltd



Veterinary Supplies



As we are the supplier to many practices for their orthopaedic and instrument needs, we thought it would be useful to add a range of consumable items linked in to orthopaedics. We are particularly excited about our collaborations with Biomendex and becoming there exclusive uk supplier of the Adaptos[®] range and also with B Braun with there Suture lines.



www.n2-uk.com



+44 (0)23 9323 3265



www.facebook.com/n2ukltd



sales@n2-uk.com



office@n2-uk.com for orders

Bone Wax

SMI Bone Wax is indicated to control bleeding on bone surfaces. It acts as a mechanical barrier to stop bleeding, exerting a haemostatic function. Composition mixture of whitened bee wax (80%) and isopropyl palmitate (20%). It is slightly yellowish and has a characteristic bee wax odour.

Will help control the bleeding form cut cancellous bone by blocking the blood vessels



Order Code	Description	RRP
BONE-WAX-12	Bone Wax - Pack 12	£32.00



Skin stapler and Remover

Indications SMI Skin Stapler can be used for a variety of skin closures

Benefits Ergonomic and user friendly design

Clear view of operating site

Easy to check remaining staples



Order Code	Description	RRP
SS-REM	Skin Stapler and Remover	£8.00

Skin Marker

Indications SMI Skin Stapler can be used for a variety of skin closures

Benefits Ergonomic and user friendly design

Clear view of operating site

Easy to check remaining staples



Order Code	Description	RRP
SKIN-MARK-6	Sterile Skin Marker with Ruler -pack 6	£15.00
SKIN-MARK-36	Sterile Skin Marker with Ruler -pack 36	£72.00



Adaptos[®] bone graft substitute

Groundbreaking tissue regeneration technology

The unique composite material of Adaptos[®] granules is a combination of the best structural and functional features of both gel-like and ceramic bone graft substitutes.

Main features

Patented 3D composite structure incorporates the osteoconductive and osteostimulative ceramic materials in a durable, elastic, and resilient carrier material resulting in superior handling qualities and operational performance.

Highly bioactive

- Bonds directly to bone tissue
- Osteostimulative bioactive glass induces bone formation
- Faster and safer healing

Press-fit implantation


- Recovers after compression and locks to defect site
- Excellent tissue contact on site
- Faster osseointegration

Resilient


- Durable and elastic composite material
- Will not slip, break or crumble even in demanding use
- Squeezable and can be stretched

Optimal porosity

- Pore structure similar to cancellous bone
- Ideal pore size for bone formation
- Allows vascular ingrowth



200µm



1mm

Other Adaptos[®] features

Adaptos[®] delivers the largely used and safe raw materials in a unique composite form: the osteostimulative and osteoconductive ceramic granules are homogenously embedded into a highly porous and resilient polymer carrier.

- Bonds directly to bone tissue
- Bioactive glass as an osteostimulative component
- Completely biodegradable
- Allows vascularization
- Provides natural structural support
- Intuitive usability



🌐 www.n2-uk.com

☎ +44 (0)23 9323 3265

🐦 @n2ukcom

📘 www.facebook.com/n2ukltd

✉ sales@n2-uk.com

🏠 office@n2-uk.com for orders



Adaptos®Vet by Biomendex

The first truly adaptable bone graft substitute for veterinary use with superior handling qualities and performance.

Order Code	Description	RRP
BG-UF-05-T	AdaptosVet Ultrafine (0.2-0.5mm)Granules 0.5cc Tube	£68.00
BG-UF-05-PK5-T	AdaptosVet Ultrafine (0.2-0.5mm)Granules 0.5cc - Tube (Pack 5)	£278.00
BG-F-10-PK5-T	AdaptosVet fine (0.5-1.0mm)Granules 1cc - Pack 5 Tube	£339.00
BG-F-05-PK5-T	AdaptosVet fine (0.5-1.0mm)Granules 0.5cc - Tube (Pack 5)	£278.00
BG-F-10-S	AdaptosVet fine (0.5-1.0mm)Granules 1cc Syringe	£96.00
BG-F-05-S	AdaptosVet fine (0.5-1.0mm)Granules 0.5cc Syringe	£79.00
BG-F-10-T	AdaptosVet fine (0.5-1.0mm)Granules 1cc Tube	£83.00
BG-F-05-T	AdaptosVet fine (0.5-1.0mm)Granules 0.5cc Tube	£68.00
BG-M-10-PK5-T	AdaptosVet Medium (1.0-2.0mm)Granules 1.0cc - Tube (Pack 5)	£339.00
BG-M-05-PK5-T	AdaptosVet Medium (1.0-2.0mm)Granules 0.5cc - Tube (Pack 5)	£278.00
BG-M-10-S	AdaptosVet Medium (1.0-2.0mm)Granules 1.0cc Syringe	£96.00
BG-M-05-S	AdaptosVet Medium (1.0-2.0mm)Granules 0.5cc Syringe	£79.00
BG-M-10-T	AdaptosVet Medium (1.0-2.0mm)Granules 1.0cc Tube	£83.00
BG-M-05-T	AdaptosVet Medium (1.0-2.0mm)Granules 0.5cc Tube	£68.00

Adaptos®Vet Ultra Fine



Recommended for



Oral & Maxillofacial

Adaptos®Vet Ultra Fine (0.2-0.5 mm)

0.5 cc options:

5x0.5 cc tube

1x0.5 cc tube

Ultra Fine granules are especially suitable for void filling of bone cysts and tooth extraction sites and for treating small fractures, especially on small veterinary patients.

Adaptos®Vet Fine



Recommended for



Orthopaedics



Oral & Maxillofacial

Adaptos®Vet Fine (0.5-1.0 mm)

0.5 cc options:

5x0.5 cc tube

1x0.5 cc syringe

1x0.5 cc tube

1.0 cc options:

5x1.0 cc tube



1x1.0 cc syringe

1x1.0 cc tube

Fine granules are especially suitable for treating fragmented fractures, filling narrow fracture gaps and for treating general bone defects in small veterinary patients.

Adaptos®Vet applications

Adaptos®Vet products are bone graft substitutes with a unique composition and bone-like structure. Resilient Adaptos®Vet gives natural assistance for regenerating bone tissue. The elastic bone graft substitute is easy to shape and use, and it can be press-fitted into a defect for faster osseointegration. Adaptos®Vet is intended for treatment of bone defects in the whole skeletal system.

Orthopaedics 	Oral & Maxillofacial 
<ul style="list-style-type: none"> • Auto or allograft extender • Void filling and osteotomy sites • Filling and reconstruction of metaphyseal bone defects • Fracture repair • Arthrodesis 	<ul style="list-style-type: none"> • Void filling, such as tooth extraction sites • Periodontal pockets and other bone loss • Fracture repair • Cyst extraction sites and other osseous defects

Adaptos®Vet Medium



Recommended for



Orthopaedics

Adaptos®Vet Medium (1.0-2.0 mm)

0.5 cc options:

- 5x0.5 cc tube
- 1x0.5 cc syringe
- 1x0.5 cc tube

1.0 cc options:

- 5x1.0 cc tube
- 1x1.0 cc syringe
- 1x1.0 cc tube

Medium granules are especially suitable for filling larger bone defects such as bone cavities, wider fracture gaps, and for treating general bone defects in larger veterinary patients.

Case Studies

Revision of a non-union: Treatment of Italian Greyhound's fractures in radius and ulna diaphysis

A 5-year-old Italian Greyhound arrived at the clinic with complicated ulna and radius fractures. The dog had been operated on 4-5 weeks before elsewhere, but the position of the leg had started to worry his owner. During examination, instability was detected in the diaphysis of the lower third in the fracture-line, which caused the change in leg position. An X-ray showed that the fracture had not started to heal properly but it had developed into a non-union. The bone heads were reattached to each other with a locking plate and Adaptos®Vet granules were implanted into the fracture line and its immediate proximity to support bone regeneration.

Follow-ups after surgery

Healing of the non-union was monitored with x-ray post op 10 days, five weeks, ten weeks, twelve weeks, 21 weeks and after 2 years. The follow-up period was uneventful and after two years bone remodeling was visible in x-ray and even though the fracture had healed, the cortex was not yet on display.



Post-op



2 years post-op

Treatment of Bengal cat's complicated supracondylar open fracture

A 2-year-old Bengal cat suffered a complicated open fracture (class II). The fracture was in the lower part of brachium (supracondylar fracture) with moderate skin damages, and he needed a surgery.

The fracture fragments were attached with 2 locking plates to stabilize the defect area. Adaptos®Vet bone graft substitute was implanted into the inward edge and into the immediate proximity of fracture line, to support the ossification of the defect.

Follow-ups after surgery

Patient was monitored with x-ray post op two weeks, seven weeks and 14 weeks. Monitoring period was uneventful. After 22 weeks the plate and screws on the inward edge were removed, and the leg's movements were observed to be good. The plate on the outer edge was left intact.



2 weeks



22 weeks post-op

MonoPlus® The Resilient One



Absorbable suture material

MonoPlus®



When wound healing requires time

- Degradation profile designed where an extended wound support is required.
- Wide range of needle-thread combinations
- Good knot pull tensile strength retention

Indications

- Orthopaedic surgery, capsular suturing
- Colon surgery
- Abdominal wall closure



Features


Structure:	Monofilament, uncoated
Chemical composition:	Poly-p-dioxanone
Diameter:	USP 7/0 – USP 2
Colour:	violet
Absorption time:	~ 35 days
Knot tensile strength retention:	~ 35 days p. i. 50 %
Mass absorption:	~ 180 – 210 days




Needle Information

Needle Angle of Curvature

Others



S 1/8 circle
V 1/4 circle
D 3/8 circle
H 1/2 circle
F 5/8 circle



Needle Curvature	Needle Body	Needle Tip	Needle Length	Further Information
S 1/8 circle	G trapezoidal micro engraved body	T trocar tip	The numbers indicate the distance in mm between the needle tip and needle end, measured on the external curve of the needle.	f thin bodied needle
V 1/4 circle	R round body	N blunt tip		ff extra thin bodied needle
D 3/8 circle	S cutting	S cutting tip		s strong needle body
H 1/2 circle	L lancet	C short, cutting tip		ss very strong needle body
F 5/8 circle		M micro-point		v break-off needle
G straight		MP microtip/precision		b black needle
P progressive/asymptotic		m micro-needle		To take off removable needle, with a straight tug of the needle holder, the needle can be removed from the suture fiber.
J hook		mV micro vascular		
SK ski				

Example one: HRT 37s = 1/2 circle round bodied needle with trocar tip, length 37 mm, strong needle body

Order Code	Description	RRP
SUT-MONO-C0023191	Monosyn Suture Undyed 2/0 (3) 70cm DS39(M) (Pack 36)	£111.00
SUT-MONO-C2022215	Monosyn Suture Violet 3/0 (2) 70cm DS24(M) RCP (Pack 36)	£113.00
SUT-MONO-C2022216	Monosyn Violet 2/0 (3) 70cm DS24 (M) RCP (Pack 36)	£126.00
SUT-MONO-C2022204	Monosyn Violet 4/0 (1.5) 70cm DS19 (M) RCP (Pack 36)	£125.00
SUT-MONO-C2022026	Monosyn Violet 2/0 (3) 70cm HR26 (M) RCP (Pack 36)	£124.00
SUT-MONO-C2022004	Monosyn Violet 4/0 (1.5) 70cm HR17 (M) RCP (Pack 36)	£118.00
SUT-MONO-C2023636	Monosyn Undyed 2/0 (3) 70cm DSMP24 (M) RCP (Pack 36)	£155.00
SUT-MONO-C2022236	Monosyn Violet 3/0 (2) 70cm HS26 (M) RCP (Pack 36)	£79.00
SUT-MONO-C2022025	Monosyn Violet 3/0 (2) 70cm HR26 (M) RCP (Pack 36)	£122.00
SUT-MONO-C2022305	Monosyn Violet 3/0 (2) 70cm HRC26 (M) RCP (Pack 36)	£138.00
SUT-MONO-C2022118	Monosyn Violet 2/0 (3) 70cm DR26 (M) RCP (Pack 36)	£117.00
SUT-MONO-C2022247	Monosyn Violet 2/0 (3) 70cm HS37S (M) RCP (Pack 36)	£101.00
SUT-MONO-C2022005	Monosyn Violet 3/0 (2) 70cm HR26 (M) RCP (Pack 36)	£122.00
SUT-MONO-C2022355	Monosyn Violet 3/0 (2) 70cm DRC22 (M) RCP (Pack 36)	£140.00


There are many more ranges available.
Please contact us if you don't quite see what you need.

Monosyn®

The Vets Choice



Absorbable suture material



2 weeks!



What's 2 weeks got to do with Monosyn®?

Monosyn® is a mid-term synthetic absorbable suture that supports the wound a whole week longer (2 weeks in total) than other well-known suture brands.

Why is that extra week important?

For patients who need a "mid-term" monofilament suture but require wound support for longer than 7 days... Vet surgeons may opt to use a Polydioxanone suture instead which provides 28 days wound support. However, this length of wound support and longer mass absorption is not necessary.

- this option is more costly
- the handling characteristics are not desirable...



The high memory can make it difficult to handle and it can be stiff and sharp in delicate soft tissues, where **Monosyn® is much softer.**

With its 14 day wound support and excellent handling Monosyn can be a better choice for the vet and for the patient



Features


Structure:	Monofilament, uncoated
Chemical composition:	Glyconate
Diameter:	USP 7/0 – USP 2
Colour:	undyed, dyed
Absorption time:	~ 14 days
Knot tensile strength retention:	~ 14 days p. i. 50%
Mass absorption:	~ 60 – 90 days




Needle Information

Needle Angle of Curvature

Others



S 1/8 circle
V 1/4 circle
D 3/8 circle
H 1/2 circle
F 5/8 circle



Needle Curvature	Needle Body	Needle Tip	Needle Length	Further Information
S 1/8 circle	G trapezoidal micro engraved body	T trocar tip	The numbers indicate the distance in mm between the needle tip and needle end, measured on the external curve of the needle.	f thin bodied needle
V 1/4 circle	R round body	N blunt tip		ff extra thin bodied needle
D 3/8 circle	S cutting	S cutting tip		s strong needle body
H 1/2 circle	L lancet	C short, cutting tip		ss very strong needle body
F 5/8 circle		M micro-point		v break-off needle
G straight		MP microtip/precision		b black needle
P progressive/asymptotic		m micro-needle		To take off removable needle, with a straight tug of the needle holder, the needle can be removed from the suture fiber.
J hook		mV micro vascular		
SK ski				

Example one: HRT 37s = 1/2 circle round bodied needle with trocar tip, length 37 mm, strong needle body

Order Code	Description	RRP
SUT-MONO-C2024306	Monoplus Suture Violet 2/0 (3) 70cm HRC26 (M) RCP (Pack 36)	£163.00
SUT-MONO-C0024101	Monoplus Suture Violet 3/0 (3) 45cm DS24 (M) (Pack 36)	£153.00
SUT-MONO-C2024305	Monoplus Suture Violet 3/0 (3) 70cm HRC26 (M) RCP (Pack 36)	£170.00
SUT-MONO-C0024517	Monoplus Suture Violet 0 (3.5) 90cm HRT37 (M) (Pack 36)	£167.00
SUT-MONO-C0024436	Monoplus Suture Violet 2/0 (3) 90cm HR26 (M) (Pack 36)	£133.00
SUT-MONO-C0024983	Monoplus Suture Violet 4/0 (1.5) 70cm DSMP24 (M) (Pack 36)	£196.00
SUT-MONO-C0024546	Monoplus Suture Violet 1 (4) 90cm HRT40S (M) (Pack 36)	£183.00
SUT-MONO-C2024135	Monoplus Suture Violet 3/0 (2) 70cm HR22 (M) RCP (Pack 36)	£138.00
SUT-MONO-C2024134	Monoplus Suture Violet 4/0 (1.5) 70cm HR22 (M) RCP (Pack 36)	£139.00
SUT-MONO-C0024333	Monoplus Suture Violet 1 (4) 90cm HRT48 (M) (Pack 36)	£261.00
SUT-MONO-C0024988	Monoplus Suture Violet 2/0 (3) 70cm DSMP24 (M) (Pack 36)	£123.00
SUT-MONO-C0024986	Monoplus Suture Violet 3/0 (2) 70cm DSMP24 (M) (Pack 36)	£151.00
SUT-MONO-C2024306	Monoplus Suture Violet 2/0 (3) 70cm HRC26 (M) RCP (Pack 36)	£163.00
SUT-MONO-C0024101	Monoplus Suture Violet 3/0 (2) 45cm DS24 (M) (Pack 36)	£153.00

We have over 5,000 suture codes, the above are our most popular Monoplus® codes. Please let us know if you have a preference for a different suture.

There are many more ranges available.
Please contact us if you don't quite see what you need.

Surgicryl® absorbable sutures

SURGICRYL®
MONOFILAMENT
POLYDIOXANONE



Surgicryl® Monofilament
Polydioxanone
Approximately 60 days
wound support

SURGICRYL®
PGA
POLYGLYCOLIC ACID



Surgicryl® PGA
Polyglycolic acid
Approximately 30 days
wound support

SURGICRYL®
MONOFAST
POLYGLYCAPRONE 25




Surgicryl® Monofast
Polyglactaprone 25
Approximately 14 days
wound support

SURGICRYL®
910
POLYGLACTINE







Surgicryl® 910
Polyglactine 910
Approximately 30 days
wound support

SURGICRYL®
MONOFAST
POLYGLYCAPRONE 25



FOIL PACKS

Round bodied taper point						Metric	USP	Length	Order number*	Reverse cutting						Metric	USP	Length	Order number*
 1/2 Circle	HR 13 mm	0.7	6/0	45 cm	16070113	 3/8 Circle	DS 12 mm					0.7	6/0	75 cm	16071512V				
	HR 17 mm	1	5/0	45 cm	16100117							1	5/0	45 cm	16101512V				
	HR 22 mm	1.5	4/0	75 cm	16150122		DS 16 mm	1	5/0	75 cm	16101516V								
		2	3/0	75 cm	16200122			1.5	4/0	75 cm	16151516V								
	HR 26 mm	2	3/0	75 cm	16200126		DS 19 mm	1	5/0	75 cm	16101519V								
		3	2/0	75 cm	16300126			1.5	4/0	75 cm	16151519V								
	HR 36 mm	3.5	0	90 cm	16350136			1.5	4/0	75 cm	16151519V								
HR 40 mm	4	1	90 cm	16400140	2			3/0	75 cm	16201519V									
							DS 26 mm	2	3/0	75 cm	16201526V								
3								2/0	75 cm	16301526V									
							DS 30 mm	3.5	0	75 cm	16351530V								


* V - violet

Type	Monofilament
Composition	Polyglactaprone 25, a copolymer made of 75% glycolide and 25% ε-caprolactone
Coating	None
Colour	Violet / Beige – undyed
Tissue reaction	Minimal
Absorption	The hydrolytic action by which the material is broken down results in total absorption in approximately 90 to 110 days. Approximately 75% of tensile strength remain after 7 days.

Order Code	Description	RRP
SUT-SURG-MF-1610-1519	Surgicryl Monofast Polyglactaprone Suture 5/0 (1) 45cm DS19 3/8 Circle Pack 12	£35.00
SUT-SURG-MF-1615-1519	Surgicryl Monofast Polyglactaprone Suture 4/0 (1.5) 75cm DS19 3/8 Circle Pack 12	£35.00
SUT-SURG-MF-1615-1519V	Surgicryl Monofast Polyglactaprone Suture 4/0 (1.5) 75cm DS19 3/8 Circle Violet Pack 12	£35.00
SUT-SURG-MF-1620-1519	Surgicryl Monofast Polyglactaprone Suture 3/0 (2) 75cm DS19 3/8 Circle Pack 12	£35.00
SUT-SURG-MF-1630-1519	Surgicryl Monofast Polyglactaprone Suture 2/0 (3) 75cm DS19 3/8 Circle Pack 12	£35.00
SUT-SURG-MF-1620-1524	Surgicryl Monofast Polyglactaprone Suture 3/0 (2) 75cm DS24 3/8 Circle Pack 12	£35.00
SUT-SURG-MF-1630-1524	Surgicryl Monofast Polyglactaprone Suture 2/0 (3) 75cm DS24 3/8 Circle Pack 12	£35.00
SUT-SURG-MF-1620-1526V	Surgicryl Monofast Polyglactaprone Suture 3/0 (2) 75cm DS26 3/8 Circle Violet Pack 12	£35.00
SUT-SURG-MF-1630-1526V	Surgicryl Monofast Polyglactaprone Suture 2/0 (3) 75cm DS26 3/8 Circle Violet Pack 12	£35.00
SUT-SURG-MF-1635-1530V	Surgicryl Monofast Polyglactaprone Suture 0 (3.5) 75cm DS30 3/8 Circle Violet Pack 12	£35.00






There are many more ranges available.
Please contact us if you don't quite see what you need.

Type	Monofilament
Composition	Polydioxanone, a polymer made from polyester poly (p-dioxanone)
Coating	None
Colour	Violet
Tissue reaction	Minimal
Absorption	The hydrolytic action by which the material is broken down results in total absorption in approximately 180 to 210 days. Approximately 50% of tensile strength remain after 42 days.




FOIL PACKS


See our website for the complete range

Round bodied taper point		Metric	USP	Length	Order number	Reverse cutting		Metric	USP	Length	Order number
 1/2 Circle	HR 13 mm	0.7	6/0	75 cm	13070113	 3/8 Circle	DS 12 mm	0.7	6/0	45 cm	13071512
	HR 17 mm	1	5/0	75 cm	13100117		DS 19 mm	1	5/0	45 cm	13101512
		1.5	4/0	75 cm	13150117			1	5/0	45 cm	13101519
	HR 22 mm	1.5	4/0	75 cm	13150122		DS 24 mm	1.5	4/0	75 cm	13151519
		2	3/0	75 cm	13200122			2	3/0	75 cm	13201519
	HR 26 mm	2	3/0	75 cm	13200126		DS 40 mm	2	3/0	75 cm	13201524
		3	2/0	75 cm	13300126			3	2/0	75 cm	13301524
	HR 30 mm	2	3/0	75 cm	13200130		HS 36 mm	3.5	0	75 cm	13351524
		3	2/0	75 cm	13300130			4	1	90 cm	13401540
	HR 36 mm	3	2/0	75 cm	13300136		 1/2 Circle	5	2	90 cm	13501540
3.5		0	90 cm	13350136	3	2/0		75 cm	13301436		
HR 40 mm		3	2/0	90 cm	13300140	3.5		0	75 cm	13351436	
HR 48 mm	4	1	90 cm	13400148			4	1	75 cm	13401436	

TENSILE STRENGTH



ABSORPTION



Order Code	Description	RRP
SUT-SURG-PD-1315-0122	Surgicryl Monofilament Polydioxanone Suture Voilet 4/0 (1.5) 75cm HR22 1/2 Circle Pack 12	£35.00
SUT-SURG-PD-1320-0122	Surgicryl Monofilament Polydioxanone Suture Voilet 3/0 (2) 75cm HR22 1/2 Circle Pack 12	£35.00
SUT-SURG-PD-1320-0126	Surgicryl Monofilament Polydioxanone Suture Voilet 3/0 (2) 75cm HR26 1/2 Circle Pack 12	£35.00
SUT-SURG-PD-1330-0126	Surgicryl Monofilament Polydioxanone Suture Voilet 2/0 (3) 75cm HR26 1/2 Circle Pack 12	£35.00
SUT-SURG-PD-1320-1524	Surgicryl Monofilament Polydioxanone Suture Voilet 3/0 (2) 75cm DS24 3/8 Circle Pack 12	£35.00
SUT-SURG-PD-1330-1524	Surgicryl Monofilament Polydioxanone Suture Voilet 2/0 (3) 75cm DS24 3/8 Circle Pack 12	£35.00
SUT-SURG-PD-1335-1524	Surgicryl Monofilament Polydioxanone Suture Voilet 0 (3.5) 75cm DS24 3/8 Circle Pack 12	£35.00
SUT-SURG-PD-1330-1436	Surgicryl Monofilament Polydioxanone Suture Voilet 2/0 (3) 75cm HS36 1/2 Circle Pack 12	£35.00
SUT-SURG-PD-1335-1436	Surgicryl Monofilament Polydioxanone Suture Voilet 0 (3.5) 75cm HS36 1/2 Circle Pack 12	£35.00
SUT-SURG-PD-1340-1436	Surgicryl Monofilament Polydioxanone Suture Voilet 1 (4) 75cm HR26 1/2 Circle Pack 12	£35.00

There are many more ranges available.
Please contact us if you don't quite see what you need.

Absorbant Op Covers

ABSORBENT OP-COVER For quality & safety Surgical drapes material 56 g/m² non breathable 3 layers Polyolefines based film 24 g/m² Rubber based adhesive (UFD technology) 2 g/m² Hydrophilic polypropylene spunbond 30 g/m²

Colour Blue Attributes: OP cover is latex free - OP cover is soft and drape-able - Liquid distributes and absorbs quickly



Order Code	Description	RRP
OP-AB-31441545-A	31 x 44cm Absorbant Op-Cover with 1.5cm x 4.5cm Aperture - Pack 225	£110.00
OP-AB-3144-3509-A	31 x 44cm Absorbant Op-Cover with 3.5cm x 9cm Aperture - Pack 225	£110.00
OP-AB-3144-0514-A	31 x 44cm Absorbant Op-Cover with 5cm x 14cm Aperture - Pack 200	£110.00
OP-AB-4474-0720-A	44 x 74cm Absorbant Op-Cover with 7cm x 20cm Aperture - Pack 100	£95.00
OP-AB-4474-0824-A	44 x 74cm Absorbant Op-Cover with 8cm x 24cm Aperture - Pack 125	£115.00

Order Code	Description	RRP
OP-AB-3045-NA	30 X 45cm Absorbant Op-Cover - Pack 225	£80.00
OP-AB-4560-NA	45 X 60cm Absorbant Op-Cover - Pack 125	£65.00
OP-AB-6090-NA	60 X 90cm Absorbant Op-Cover - Pack 75	£58.00
OP-A-90120-NA	90 X 120cm Absorbant Op-Cover - Pack 40	£50.00
OP-120120-NA	120 X 120cm Absorbant Op-Cover - Pack 30	£50.00
OP-120180-NA	120 X 180cm Absorbant Op-Cover - Pack 15	£40.00
OP-180250-NA	180 X 250cm Absorbant Op-Cover - Pack 15	£55.00

Heamostatic Sponge

Our gelatin sponges are surgical haemostatic sponges made from first grade highly purified gelatin material. SMI-Spon is intended for use in various surgical procedures and quickly and effectively stops bleeding fast. It reduces intraoperative bleeding and surgery time.



Order Code	Description	RRP
SP-AH-805010-10	Absorbable Haemostatic Gelatin Sponge Standard - 80 x 50 x 10mm (Pack 10)	£22.00
SP-AH-80501-20	Absorbable Haemostatic Gelatin Sponge Special - 80 x 50 x 1mm (Pack 20)	£20.00
SP-AH-705010-10	Absorbable Haemostatic Gelatin Sponge Regular - 70 x 50 x 10mm (Pack 10)	£28.00
SP-AH-101010-32	Absorbable Haemostatic Gelatin Sponge Dental Cube - 10 x 10 x 10mm (Pack 32)	£15.00

Surgibond

Surgibond Vet Glue has the following uses

- provides additional suture support where required
- stops bleeding
- helps protect against infection
- good visibility due to blue colouring
- acts within seconds
- no need for sutures or dressings



Order Code	Description	RRP
VET-GLUE	Surgibond Vet Glue 2.5ml	£15.00
VET-GLUE-APP-6	Surgibond Vet Glue Applicators - pack 6	£1.50