



N2 (UK) Ltd



Veterinary Supplies

# T-Plates with Locking and Evolox<sup>®</sup>



Order Code	Description	Length	RRP
LPL-T-20-5-35	T-Plate Locking 2.0mm 2+3	35mm Long	£65.00
LPL-T-20-6-35X	T-Plate Locking 2.0mm 2+4	35mm Long	£65.00
EV-T-2427-7-48	T-Plate Evolox <sup>®</sup> 2.4/2.7mm 3+4	48mm Long	£70.00
EV-T-2427-8-55	T-Plate Evolox <sup>®</sup> 2.4/2.7mm 3+5	55mm Long	£75.00
EV-T-35-7-70	T-Plate Evolox <sup>®</sup> 3.5mm 3+4	70mm Long	£70.00
EV-T-35-8-80	T-Plate Evolox <sup>®</sup> 3.5mm 3+5	80mm Long	£80.00
EV-T-35-10-100	T-Plate Evolox <sup>®</sup> 3.5mm 3+7	80mm Long	£90.00

We are pleased to introduce the expansion of our Evolox<sup>®</sup> and Locking T Plates with emphasis of improving on current Implants.

Designed, Developed and Manufactured in the UK.

- Each Plate size can take the relevant ISO recognised Locking and Cortical Screws.
- Scalloped Edges for easier anatomical handling.
- Arthrodesis Wire hole to aid placement.
- Grooved Feature for Plate Clamping.
- Evolox<sup>®</sup> for Poly-axial Locking for better pull-out resistance and Secure Fixation.
- Can be contoured.



N2 (UK) Ltd



# T-Plates with Evolox<sup>®</sup>/Locking

## Veterinary Supplies



LPL-T-20-5-35




LPL-T-20-6-35X




EV-T-2427-7-48




EV-T-2427-8-55



EV-T-35-7-70



EV-T-35-8-80



EV-T-35-10-100

Drawings at 1:1 scale

- Made with Human Grade Stainless Steel to ISO 5832-1.
- Non-Reflective surface.
- 2mm Plate will take either 2mm Cortical, 2mm Locking or 2.2mm Cancellous.
- 2.4mm/2.7mm will take 2.4/2.7 Cortical, 2.4/2.7 Locking and 3mm Cancellous.
- 3.5mm will take either 3.5mm Cortical, 3.5mm Locking or 4.0mm Cancellous.
- Can be used with N2's Hybrid Locking Screws where Bone Stock is Limited.



- 🌐 [www.n2-uk.com](http://www.n2-uk.com)
- ☎ +44 (0)23 9323 3265
- 🐦 @n2ukcom
- 📘 [www.facebook.com/n2ukltd](http://www.facebook.com/n2ukltd)
- ✉ [sales@n2-uk.com](mailto:sales@n2-uk.com)
- 🏠 [office@n2-uk.com](mailto:office@n2-uk.com) for orders

N2(UK) Ltd is a privately Owned Independent Company.