



Hybrid Locking Kit





Veterinary Supplies

Our Hybrid Locking screws are getting more and more popular with our customers due to their usefulness in many situations and are a great product to have sterilised in your practice for those unexpected cases.

On top of this we now have the new locking adapter. This screws into a 3.5mm Locking/Evolox plate and allows you to use a 2.4mm or 2.7mm screw in any 3.5mm plates polyaxially.

It also has a slot on the washer to allow you to remove it for any reason.

The Complete kit is available preloaded with a full set of the hybrid screws in a premium box along with 2 of the Locking adapters and one of each required Locking drill Guides.



Order Code	Description	RRP
PM-RACK-HYBRID-LOCK	Empty Hybrid Premium Screw Rack	£300.00
KIT-RACK-HYBRID-LOCK	Full Hybrid Screw Washer kit in Premium Screw Rack	£2,000.00
LOC-AD-35-2427	Locking Adapter -3.5mm down to 2.4/2.7mm	£20.00



- 🌐 www.n2-uk.com
- ☎ +44 (0)23 9323 3265
- 🐦 @n2ukcom
- 📘 www.facebook.com/n2ukltd
- ✉ sales@n2-uk.com
- 📧 office@n2-uk.com for orders